

CPAP MACHINE	DISPOSABLE FILTER
EVERY 5 YEARS	2 TIMES A MONTH
NASAL CUSHION/PILLOWS	FULL FACE CUSHION
2 TIMES A MONTH	EVERY 30 DAYS
FRAME	TUBING
EVERY 3 MONTHS	EVERY 3 MONTHS
HEADGEAR	WATER CHAMBER

Keep yourself healthy and infection free by replacing your CPAP supplies on schedule

919.462.8081 | 130 Preston Executive Drive Cary, NC | parkwaysleep.com